

# TECHNICAL DATA SHEET

## « NG FO-MT+ »

Tannins for Wines, Beers and Spirits

#### **Composition:**

The Next Generation FO-MT+ is a Tannin extracted from French Oak that has been Open-Air seasoned for a minimum of 24 Months and designed for use with both Wines and Spirits. This product is available in Liquid and Freeze-Dried formulations.

#### **General Characteristics:**

The NG FO-MT+ is an ellagic tannin extracted with our properitary aqueous method from both Quercus Robur and Quercus Petraea. Next Generation FO-MT+ is applicable for White, Red and Rose Wines during ageing as well as many different types of Spirits.

Technological assistant - Oenological product / Customs code = 32019025 - Tannin.

### **Applications:**

Ageing:

- Bring new barrel profile French Oak / Medium Toast +
- Increase sucrosity & length on the palate
- Add tannins

#### Before bottling:

- Bring aged character
- Adjust the profile of the wine with a fast integration

#### **Dosage Recommendation:**

- White Wines: 50 - 400mL/hL

- Red Wines & rosés : 100 - 500mL/hL

- Spirits range: 500 - 2000mL/hL

These figure are indicative, a laboratory test must be carried out, in order to validate the dose of use.

#### **Compliance:**

Product approved for Distilling and Winemaking by the TTB.

Legal Limit: The residual amount of tannin shall not exceed 0.8 g/L in White Wine and 3.0 g/L in Red Wine (in gallic acid). Only tannin which does not impart color may be used.

Total tannin shall not be increased by more than 150 mg/L (in tannic acid).

For Spirits, please consult legal limits for each individual application.

#### **Efficacy:**

Liquid: 6 Months under refrigeration between 35 and 41° Farenheit.

Freeze Dried: 3 Years in Cellar conditions.

## Packaging Available:


Liquid: 1L +

Freeze Dried: 1/4 kg +

# TECHNICAL DATA SHEET

## « NG FO-MT+ »

## **Spider Graph:**



This product is proudly represented in North America by:



www.cheneetfils.com

Contact:

Dan Sacks: 707-933-6477 <u>dan@cheneetfils.com</u>